

Composition:

Each 1 ml of salinose PLUS contains:

3.3 mg of Sodium Chloride

8.3 mg of Sodium Bicarbonate

Indications:

Isotonic saline solution to relieve nasal congestion.

Dosage and Administration:

Dosage:

- Adults and children (from 2 years onwards): spray 1-2 puffs into each nostril several times daily.

Application and Administration:

1- Push down the white actuator one time or several times until the first puff is dispensed. Now the spray is ready to use.

2- Insert the spray nozzle into the nostril and apply pump mechanism. In doing so keep the bottle vertical and don't bend head backwards. During the spraying process breathe easily through the nose.

Contraindications:

There are no experience of contraindications against salinose PLUS, however after nose operations or accidents a doctor should be consulted before using the product.

Side Effects:

If used correctly, side effects are not to be expected.

Warnings and Precautions:

- Due to reasons of hygiene and to avoid infections, every bottle of salinose PLUS should only be used by one person.
- After each application, the spray nozzle should be wiped and cleaned and the protective cap should be replaced on the bottle again.
- For nasal use only.
- For pregnant and nursing mothers, please consult your physician before using salinose PLUS.

Storage Conditions:

- Store below 30° C.

Presentation:

salinose PLUS is available in 75 ml jet spray.

Produced By :

Middle East Pharmaceutical Industries Co. Ltd (Avalon Pharma), Riyadh, Saudi Arabia.
Tel +966(11)2653948

Revision Date:01/2023